

Invisible below
the temporary
restoration

Visible during
removal of
excess
material

Bifix® Temp

DUAL-CURING, COMPOSITE-BASED TEMPORARY
LUTING MATERIAL



Bifix® Temp

EASE OF USE AND OUTSTANDING AESTHETICS COMBINED IN A SINGLE PRODUCT

With Bifix Temp the perfect aesthetic result is already achieved when luting temporary restorations. The translucent, tooth-like Universal shade is ideal for highly aesthetic or very thin and thus transparent temporary restorations, does not shine through them and therefore enhances the natural appearance of the restoration. What's more, it is also extremely easy to use.

Indication


Non-permanent luting of temporary restorations

Advantages


- Film thicknesses of less than 5 µm – the accuracy of fit is not affected
- Dual-curing – flexible setting times by use of a polymerisation light
- Simple removal of excess material – can be removed in one piece following brief light-curing
- In the tried-and-tested QuickMix syringe – automatic mixing without errors
- Short mixing tip – minimal material loss



Influence of the luting material on Structur 3



Opaque temporary luting materials can have a major effect on the aesthetics of temporary restorations being fitted. In the anterior region in particular this is a highly undesirable factor.



Results of survey among testing dentists

REF 1055 QuickMix syringe 5 ml universal, mixing tips type 9

REF 2188 Mixing tips type 9, 50 pcs.

VOCO GmbH
 Anton-Flettner-Straße 1-3
 27472 Cuxhaven
 Germany
 Tel.: +49 (0) 4721-719-0
 Fax: +49 (0) 4721-719-140
 info@voco.com
 www.voco.com

Available from:

